



Linking people to care

## 9. Support Documents and Forms



**Support Documents and Forms**

Checklist for Conducting Rapid HCV Tests ..... 9-2

HCV Positive Enhance Risk Assessment Tool ..... 9-4

Client Test Result Log ..... 9-5

Product Information Training ..... 9-7

Proficiency Testing Panel Results ..... 9-12

Investigational and Remedial Action on Unacceptable Proficiency Testing ..... 9-13

Clinic Communication and Complaint Log ..... 9-14

# Checklist for Conducting Rapid HCV Tests

---

## General Test Preparation

- Conditions for testing verified (temperature and lighting)
- Clock or timer made available
- Expiration date verified on Pouch
- Clean, disposable absorbent workspace cover
- Manufacturer's Stand used
- Stand is on flat level surface
- Test Device left in pouch until needed (not contaminated)
- Absorbent Packet included in Pouch
- Vial slid into Stand
- Vial is completely seated in Stand
- Two holes in back of Test Device not covered

## Sample Collection and Test Procedures - Fingerstick Whole Blood Collection

- Disposable gloves worn
- Client finger cleaned with antiseptic wipe
- Finger allowed to dry thoroughly
- Second drop of blood collected
- Loop was completely filled with blood
- Loop stirred into Developer Solution
- Solution turned pink
- Pad on the Test Device touched the bottom of the vial
- Results Window faced forward
- Timer used
- Results read between 20 and 40 minutes after Test Device inserted into the vial
- Test Device read while in Developer Solution Vial
- Test results properly recorded

# Checklist for Conducting Rapid HCV Tests

---

## Sample Collection and Test Procedures - Venipuncture Whole Blood Collection

- Specimen collection conducted with an EDTA, lithium heparin, sodium heparin or sodium citrate test tube
- Loop was completely filled with blood
- Loop stirred into Developer Solution
- Solution turned pink
- Pad on the Testing Device touched the bottom of the vial
- Results Window faced forward
- Timer used
- Results read between 20 and 40 minutes after Testing Device inserted into the vial
- Test Device read while in Developer Solution Vial

## Quality Control (To be conducted as per manufacturer's guidelines)

- Kit Controls run; Date and Time\_\_\_\_\_
- Kit Controls verify Control Test Results match the expected results
- Good lighting used with Controls
- Test viewed good lighting conditions
- Clinic room temperature checked and recorded

**Comments:** \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

# HCV Positive Enhanced Risk Assessment Tool

---

1. Did you receive a blood transfusion or solid organ transplant (heart, lung, liver, pancreas, kidney) before July 1992?  No  Yes
2. Did you receive clotting factor concentrates produced before 1987?  No  Yes
3. Have you ever received hemodialysis?  No  Yes
4. Have you had blood tests that showed a liver problem?  No  Yes
5. Have you had a needlestick injury working in a health care setting?  No  Yes
6. Did your mother have hepatitis C when you were born?  No  Yes
7. Have you shared a toothbrush, razor, or any other item that might have blood on it (visible or not) with a person who has hepatitis C?  No  Yes
8. Have you or any of your sex partner(s) injected illegal drugs, even if it was only one time many years ago?  No  Yes
9. Have you ever been told by a medical provider that you have HIV infection?  No  Yes

Clinic Name/Site Location: \_\_\_\_\_

Name/ID Number: \_\_\_\_\_ Date: \_\_\_\_\_

Counselor/Tester: \_\_\_\_\_ Date: \_\_\_\_\_

Comments: \_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_

## OraQuick® HCV Rapid Antibody Test Client Test Result Log

Clinic Name: \_\_\_\_\_








Testing Location: \_\_\_\_\_

Client ID/Name	Counselor Code/Initials	Test Dates mm/dd/yy	Test # Performed	Lot #/Exp. Date of Test	Time Test Performed	Testing Room Temperature	Time Test Interpreted	Testing Room Temperature	Test Result	Date Client Received Result	Result & Type of Supplemental Serum or Plasma	Date Client Received Test/Suppl Result	How Result Received by: Phone (PH) or Result in Person (IP)
			<input type="checkbox"/> First <input type="checkbox"/> Second		____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	 C T Non-Reactive Reactive (circle one)		<input type="checkbox"/> Serum <input type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative		PH IP (circle one)
			<input type="checkbox"/> First <input type="checkbox"/> Second		____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	 C T Non-Reactive Reactive (circle one)		<input type="checkbox"/> Serum <input type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative		PH IP (circle one)
			<input type="checkbox"/> First <input type="checkbox"/> Second		____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	 C T Non-Reactive Reactive (circle one)		<input type="checkbox"/> Serum <input type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative		PH IP (circle one)
			<input type="checkbox"/> First <input type="checkbox"/> Second		____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	 C T Non-Reactive Reactive (circle one)		<input type="checkbox"/> Serum <input type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative		PH IP (circle one)
			<input type="checkbox"/> First <input type="checkbox"/> Second		____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	 C T Non-Reactive Reactive (circle one)		<input type="checkbox"/> Serum <input type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative		PH IP (circle one)
			<input type="checkbox"/> First <input type="checkbox"/> Second		____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	____ a.m. ____ p.m. (circle one)	____ °C ____ °F (circle one)	 C T Non-Reactive Reactive (circle one)		<input type="checkbox"/> Serum <input type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative		PH IP (circle one)

## OraQuick® HCV Rapid Antibody Test Client Test Result Log

Clinic Name: \_\_\_\_\_

Testing Location: \_\_\_\_\_

Client ID/Name	Counselor Code/Initials	Test Dates mm/dd/yy	Test # Performed	Lot #/Exp. Date of Test	Time Test Performed	Testing Room Temperature	Time Test Interpreted	Testing Room Temperature	Test Result	Date Client Received Result	Result & Type of Supplemental Serum or Plasma	Date Client Received Test/Suppl Result	How Result Received by: Phone (PH) or Result in Person (IP)
Virginia Hall	#938-LB	09/05/11	<input checked="" type="checkbox"/> First <input type="checkbox"/> Second	#1234567	8:45 a.m. (circle one)	72.6 °F (circle one)	9:08 a.m. (circle one)	72.6 °F (circle one)	 Reactive (circle one)	9/5/11	<input checked="" type="checkbox"/> Serum <input checked="" type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative	N/A	PH IP (circle one)
#12345678	#938-LB	09/05/11	<input checked="" type="checkbox"/> First <input type="checkbox"/> Second	#1234567	9:42 a.m. (circle one)	72.4 °F (circle one)	10:10 a.m. (circle one)	72.0 °F (circle one)	 Reactive (circle one)	Positive 9/5/11	<input type="checkbox"/> Serum <input checked="" type="checkbox"/> Plasma <input checked="" type="checkbox"/> Positive <input type="checkbox"/> Negative	9/16/03	PH IP (circle one)
Mark Dewitt	#938-LB	09/05/11	<input checked="" type="checkbox"/> First <input type="checkbox"/> Second	#1234567	11:55 a.m. (circle one)	72.0 °F (circle one)	12:17 a.m. (circle one)	72.0 °F (circle one)	 Reactive (circle one)	Invalid N/A	<input type="checkbox"/> Serum <input type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative		PH IP (circle one)
Theresa Howard	#938-LB	09/05/11	<input checked="" type="checkbox"/> First <input type="checkbox"/> Second	#1234567	12:30 a.m. (circle one)	72.1 °F (circle one)	12:50 a.m. (circle one)	72.2 °F (circle one)	 Reactive (circle one)	9/5/11	<input checked="" type="checkbox"/> Serum <input checked="" type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative	N/A	PH IP (circle one)
#12345680	#938-LB	09/05/11	<input type="checkbox"/> First <input checked="" type="checkbox"/> Second	#1234567	12:30 a.m. (circle one)	72.1 °F (circle one)	12:50 a.m. (circle one)	72.2 °F (circle one)	 Reactive (circle one)	9/5/11	<input type="checkbox"/> Serum <input type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative	N/A	PH IP (circle one)
Theresa Howard	#938-LB	09/05/11	<input type="checkbox"/> First <input checked="" type="checkbox"/> Second	#1234567	12:30 a.m. (circle one)	72.1 °F (circle one)	12:50 a.m. (circle one)	72.2 °F (circle one)	 Reactive (circle one)	9/5/11	<input type="checkbox"/> Serum <input type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative	N/A	PH IP (circle one)
#12345680	#938-LB	09/05/11	<input type="checkbox"/> First <input checked="" type="checkbox"/> Second	#1234567	12:30 a.m. (circle one)	72.1 °F (circle one)	12:50 a.m. (circle one)	72.2 °F (circle one)	 Reactive (circle one)	9/5/11	<input type="checkbox"/> Serum <input type="checkbox"/> Plasma <input type="checkbox"/> Positive <input type="checkbox"/> Negative	N/A	PH IP (circle one)

# OraQuick® HCV Product Information Training

---

Clinic Name/Site Location: \_\_\_\_\_

Name: \_\_\_\_\_ Date: \_\_\_\_\_

Score: \_\_\_\_\_ Trainer/Tester: \_\_\_\_\_

Select the best response for statements 1 through 6:

- 1. The complete storage temperature range of the OraQuick® HCV Test Kit is**
  - (a) 2–8°C; 36–46°F
  - (b) 2–30°C; 36–86°F
  - (c) comfortable room temperature - 15°–37°C; 59°–99°F
  
- 2. The complete storage temperature range of the OraQuick® HCV Visual Reference Panel is**
  - (a) 15–30°C; 59–86°F
  - (b) 2–30°C; 36–86°F
  - (c) comfortable room temperature - 15°–37°C; 59°–99°F
  
- 3. The complete storage temperature range of the OraQuick® HCV Kit Controls is**
  - (a) 2–8°C; 36–46°F
  - (b) 2–30°C; 36–86°F
  - (c) comfortable room temperature - 15°–37°C; 59°–99°F
  
- 4. The acceptable temperature range for performing OraQuick® HCV Test is**
  - (a) 2–8°C; 36–46°F
  - (b) 2–37°C; 35–99°F
  - (c) comfortable room temperature - 15°–37°C; 59°–99°F
  
- 5. According to the manufacturer's instructions, the acceptable time to read the OraQuick® HCV Test Device result is**
  - (a) 10 to 30 minutes
  - (b) 20 to 40 minutes
  - (c) 20 to 80 minutes
  - (d) 10 to 60 minutes
  
- 6. The three possible OraQuick® HCV Test Device result outcomes are**
  - (a) reactive, non-reactive, borderline
  - (b) reactive, non-reactive, inconclusive
  - (c) reactive, non-reactive, weakly reactive
  - (d) reactive, non-reactive, invalid
  
- 7. The blood-filled Specimen Collection Loop**
  - (a) should be rapidly dipped in the Developer Solution Vial and discarded
  - (b) should be stirred in the Developer Solution Vial and then discarded
  - (c) can be left in the Developer Solution Vial for up to 10 minutes and discarded
  - (d) should be stirred in the Developer Solution Vial and then saved until the test is complete



Select True or False for Statements 7 – 16:

**8. When conducting an OraQuick® HCV Test Control, if the positive and/or negative control does not give the correct result(s), clients can still be tested with the OraQuick® HCV Kits.**

True                      False

**9. If the absorbent packet is not present when opening the OraQuick® HCV pouch, the pouch contents should be allowed to remain open for 5 – 10 minutes before using.**

True                      False

**10. The 2 holes in the back of the OraQuick® HCV Test Device must be covered after placing the device into the Developer Solution Vial.**

True                      False

**11. The built-in procedural control in the OraQuick® HCV Test Device is intended to confirm that the patient sample has moved past the Test (T) area.**

True                      False

**12. The Developer Solution Vial must turn pink after adding the fingerstick whole blood sample.**

True                      False

**13. The OraQuick® HCV Test Device should not be removed from the Developer Solution Vial before reading the Test result.**

True                      False

**14. The first drop of blood from a fingerstick can be use to perform the OraQuick® HCV test.**

True                      False

**15. An OraQuick® HCV reactive test result is interpreted as positive test for the presence of HCV antibodies.**

True                      False

**16. OraQuick® HCV is currently CLIA-waiver approved for use in the U.S. with fingerstick whole blood specimens *only*.**

True                      False

# OraQuick® HCV Product Information Training Answer Key

---

- 1) **B.** Storage conditions for the OraQuick® HCV Test Kits is 2–30°C; 36–86°F.
- 2) **A.** Storage conditions for the OraQuick® HCV Visual Reference Panel is 15–30°C; 59–86°F.
- 3) **A.** Storage conditions for the OraQuick® HCV Kit Controls is 2–8°C; 35–46°F.
- 4) **C.** Acceptable temperature range for performing the OraQuick® HCV Test Device is a comfortable room temperature of 15°–37°C; 59°–99°F.
- 5) **B.** Acceptable times to read the OraQuick® HCV Test Device is 20 minutes to 40 minutes. **DO NOT** attempt to read the test result after the 40 minute test development.
- 6) **D.** The three possible OraQuick® HCV Test result outcomes is (R) Reactive, (NR) Non-Reactive and (INV) Invalid. An Invalid test result cannot be interpreted and is an indication that there was a problem running the test, either related to the specimen or the Device. A repeat test should be performed with a new Pouch and new sample.
- 7) **B.** The blood-filled Specimen Collection Loop should be gently stirred in the Developer Solution Vial and then immediately discarded in a bio-hazard waste container.
- 8) **False.** If the kit controls DO NOT produce the expected outcome of a reactive and non-reactive result prior to testing a patient sample, re-run the controls. If the test devices still do not produce a reactive and non-reactive result, contact OraSure Technologies. DO NOT conduct any patient testing.
- 9) **False.** An absorbent packet is included in each Device Pouch to ensure that moisture levels are maintained and do not compromise the Test Device performance. The absence of an absorbent packet means that the Test Device may have been compromised during storage. Immediately discard package, device and developer vial and open a new pouch to proceed.
- 10) **False.** The two holes at the back of the OraQuick® HCV Test Device are part of the design of the lateral flow system. Blocking these holes with labels and or other materials will interfere with the test development.
- 11) **True.** The built-in procedural control on the OraQuick® Test Device is designed to verify that the chemistry of the test has flowed past the “T” or Test Line of the Device and that a sample was added to the Developer Vial Solution.
- 12) **True.** If the blood-filled loop has been properly introduced to the Developer Solution Vial, the solution will turn a shade of pink indicating sample has been mixed.
- 13) **True.** The OraQuick® HCV Test Device should not be removed from the Developer Solution Vial until the test result has been read, interpreted and documented.
- 14) **False.** The first droplet of blood should be wiped away from the finger. Typically, the first droplet will contain tissue sample as well as blood. This may interfere with the test performance. Apply the loop to the second droplet for a clean sample.
- 15) **True.** An OraQuick® HCV Reactive test results is reported as a POSITIVE. OraQuick® HCV Rapid Antibody Test is a qualitative immunoassay test. A positive test result indicates that the presence of HCV antibodies have been detected. The patient should undergo appropriate clinical follow-up for supplemental testing.
- 16) **False.** The OraQuick® HCV Rapid Antibody Test is CLIA-waived approved in fingerstick **AND** venipuncture whole blood.

# OraQuick® HCV Rapid Test Result Panel Training

Clinic Name/Site Location: \_\_\_\_\_ Date: \_\_\_\_\_

Name: \_\_\_\_\_ Trainer/Tester: \_\_\_\_\_

Score: \_\_\_\_\_

1



Result: \_\_\_\_\_

2



Result: \_\_\_\_\_

3



Result: \_\_\_\_\_

4



Result: \_\_\_\_\_

5



Result: \_\_\_\_\_

6



Result: \_\_\_\_\_

7



Result: \_\_\_\_\_

8



Result: \_\_\_\_\_

9



Result: \_\_\_\_\_

10



Result: \_\_\_\_\_

Write the Result on the line below each Test Device: Non-Reactive (NR); Reactive (R); Invalid (INV)

# OraQuick® HCV Test Result Panel Training Answer Key

---

## OraQuick® Rapid Test Result Panel Training - Answer Key

- 1) **NR.** Test line appears only in the “C” designated area, indicating no detection of HCV antibodies are present at the time the test was conducted. Patient is presumed not to be infected with HCV.
- 2) **R.** Test lines appear in the “C” and “T” designated areas, indicating that the presence of HCV antibodies have been detected in the specimen. The patient is presumed to be infected with HCV. The strength of the “T” line has no direct correlation to a quantitative interpretation of the HCV virus. A reactive result in the OraQuick® HCV Rapid Antibody Test should undergo appropriate clinical follow-up, according to CDC recommendations for supplemental testing.
- 3) **INV.** No test lines appear in the “T” or “C” designated areas. No interpretation can be made. Repeat the test with a new Pouch and a new specimen.
- 4) **R.** Test lines appear in the “C” and “T” designated areas, indicating that the presence of HCV antibodies have been detected in the specimen. The patient is presumed to be infected with HCV. The strength of the “T” line has no direct correlation to a quantitative interpretation of the HCV virus. A reactive result in the OraQuick® HCV Rapid Antibody Test should undergo appropriate clinical follow-up, according to CDC recommendations for supplemental testing.
- 5) **INV.** No test lines appear in the “C” designated area and only a partial line appears in the “T” designated area. This result means there was a problem running the test, either related to the specimen or to the Device. No interpretation can be made. Repeat the test with a new Pouch and a new sample.
- 6) **R.** Test lines appear in the “C” and “T” designated areas, indicating that the presence of HCV antibodies have been detected in the specimen. The patient is presumed to be infected with HCV. The strength of the “T” line has no direct correlation to a quantitative interpretation of the HCV virus. A reactive result in the OraQuick® HCV Rapid Antibody Test should undergo appropriate clinical follow-up, according to CDC recommendations for supplemental testing.
- 7) **INV.** No test lines appear in the “T” or “C” designated areas. The test result window has not cleared revealing the test result rendering it impossible to read or interpret. This result means there was a problem running the test, either related to the specimen or to the Device. No interpretation can be made. Repeat the test with a new Pouch and a new sample.
- 8) **INV.** Partial test lines appear inside the “T” and “C” designated areas. While on appearance, it would seem that a POSITIVE could be interpreted, the partial development of the test lines are not fully developed within the designated areas. No interpretation can be made. Repeat the test with a new Pouch and a new sample.
- 9) **NR.** Test line appears only in the “C” designated area, indicating no detection of HCV antibodies are present at the time the test was conducted. Patient is presumed not to be infected with HCV.
- 10) **R.** Test lines appear in the “C” and “T” designated areas, indicating that the presence of HCV antibodies have been detected in the specimen. The patient is presumed to be infected with HCV. The strength of the “T” line has no direct correlation to a quantitative interpretation of the HCV virus. A reactive result in the OraQuick® HCV Rapid Antibody Test should undergo appropriate clinical follow-up, according to CDC recommendations for supplemental testing.

# Proficiency Testing Panel Results

---

## Specimen Proficiency Panel Interpretation:

For each specimen, indicate the result with a checkmark for either REACTIVE, NON-REACTIVE, or INVALID.

Assay Lot #: \_\_\_\_\_

Specimen A:  REACTIVE  NON-REACTIVE  INVALID

Specimen B:  REACTIVE  NON-REACTIVE  INVALID

Specimen C:  REACTIVE  NON-REACTIVE  INVALID

Specimen D:  REACTIVE  NON-REACTIVE  INVALID

Specimen E:  REACTIVE  NON-REACTIVE  INVALID

To be completed by study monitor:

<input type="checkbox"/> CORRECT	<input type="checkbox"/> INCORRECT
<input type="checkbox"/> CORRECT	<input type="checkbox"/> INCORRECT
<input type="checkbox"/> CORRECT	<input type="checkbox"/> INCORRECT
<input type="checkbox"/> CORRECT	<input type="checkbox"/> INCORRECT
<input type="checkbox"/> CORRECT	<input type="checkbox"/> INCORRECT

Clinic Name/Site Location: \_\_\_\_\_

Name: \_\_\_\_\_ Date: \_\_\_\_\_

Score: \_\_\_\_\_ Trainer/Tester: \_\_\_\_\_

Comments: \_\_\_\_\_  
\_\_\_\_\_  
\_\_\_\_\_

# Investigational and Remedial Action on Unacceptable Proficiency Testing

---

Date of Investigation: \_\_\_\_\_ Clinic Name/Site Location: \_\_\_\_\_

Prepared by: \_\_\_\_\_

Client Sample: <input type="checkbox"/> Yes <input type="checkbox"/> No	Control Sample: <input type="checkbox"/> Yes <input type="checkbox"/> No
Date of Testing:	Time of Testing:
Lot # of Test Device:	Exp. Date of Test Device:
Unacceptable (Reported) Result:	
Acceptable Result Range:	
Day of Testing - Quality Control Results Reviewed: <input type="checkbox"/> Yes <input type="checkbox"/> Acceptable <input type="checkbox"/> Not Acceptable Indicate Corrective Action:	
Clerical/Transcription Review: <input type="checkbox"/> Acceptable <input type="checkbox"/> Not Acceptable Indicate Corrective Action:	
Was Patient Reported Results Affected? <input type="checkbox"/> No, ( <i>skip to next section</i> ) <input type="checkbox"/> Yes Indicate Corrective Action:	
Classification of Problem: <input type="checkbox"/> Clerical <input type="checkbox"/> Technical <input type="checkbox"/> Methodology <input type="checkbox"/> Problem with Client  <input type="checkbox"/> Training Issue <input type="checkbox"/> No Explanation	
Conclusions:	
Corrective Actions/System Change(s) To Prevent Recurrence:	
Supervisor:	Date:
Lab Director:	Date:

Upon Completion - This Record Must be Maintained According to Local Regulations

## Clinic Communication and Complaint Log

Clinic Name/Site Location:		
Date Reported:	Time:	Initiated By:
Source of Communication/Complaint:		
Date of Occurrence:	Time:	
Narrative of Event (If necessary):		
Immediate Corrective Action Taken:		
Does the written procedure cover how to deal with this event? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Not Applicable		
<b>*If No – Procedure must be updated within fifteen days from date of event.</b>		
<b>If Yes – Was the written procedure followed?</b> <input type="checkbox"/> Yes <input type="checkbox"/> No	<b>If No – Why not? Explain Below</b>	
<b>Follow-up Activities Required?</b> <input type="checkbox"/> Yes <input type="checkbox"/> No	<b>If Yes-Indicate what and date to be completed below</b>	
Form Completed by:		Date:
Signature		Date:

Upon Completion - This Record Must be Maintained According to Local Regulations

## Clinic Communication and Complaint Log

Clinic Name/Site Location: <b>Schnectady Women's Health Clinic - Schnectady, New York</b>		
Date Reported: <b>9/5/11</b>	Time: <b>9:50 a.m.</b>	Initiated By: <b>Josephine Parker</b>
Source of Communication/Complaint: <b>Kathy DeWitt's (#1234581) Rapid HCV Screening Test revealed a "Positive" test result. Client verbally provided Informed Consent but refuses clinical follow-up and supplemental testing.</b>		
Date of Occurrence: <b>9/5/11</b>	Time: <b>9:00 a.m.</b>	
Narrative of Event (If necessary): <b>Explained to Ms. DeWitt's the importance of clinical follow-up appointments and the need to perform supplemental testing for the detection of viremia as well as including importance of receiving future medical care and treatment. Reviewed availability of local programs and further counseling information. Re-emphasized the issues of protection from potential exposure to partner.</b>		
Immediate Corrective Action Taken:		
Does the written procedure cover how to deal with this event? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> Not Applicable <b>*If No – Procedure must be updated within fifteen days from date of event.</b>		
If Yes – Was the written procedure followed? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No	If No – Why not? Explain Below	
Follow-up Activities Required? <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No	If Yes-Indicate what and date to be completed below	
Form Completed by: <b>Josephine Parker</b>		Date: <b>9/5/11</b>
Signature <i>Diane Rancer - Supervisor</i>		Date: <b>9/5/11</b>

Upon Completion - This Record Must be Maintained According to Local Regulations





**OraSure Technologies**

220 East First Street  
Bethlehem, PA 18015  
1-800-ORASURE (672-7873)  
[www.OraSure.com](http://www.OraSure.com)

©2011, 2012 OraSure Technologies, Inc.  
Item#HCV0068 (rev. 04/12)